



Protocol

In-Gel Digestion by Trypsin

- Excise spots/bands cutting as close to the edge of the spot/band as possible (to reduce the volume of 'background' gel)
- Transfer gel plug into a reaction tube (Eppendorf)
- Prepare the following solutions:

25 mM NH₄HCO₃ (100 mg/50 ml)

25 mM NH₄HCO₃ in 50% ACN

10 mM DTT in 25 mM NH₄HCO₃ (1.5 mg/mL)

55 mM iodoacetamide in 25 mM NH₄ HCO₃ (10 mg/mL)

50% ACN/5% formic acid (may substitute TFA or acetic acid)

12.5 ng/μL trypsin in 25mM NH₄HCO₃ (freshly diluted)

- Add ~100μL of 25mM NH₄HCO₃/50% ACN for 10 min
- Remove the supernatant and discard.
- Repeat steps 4 and 5 once or twice.
- Add 25 μL 10 mM DTT in 25 mM NH₄HCO₃. Allow reaction to proceed at 57 °C for 1 hr.
- Remove supernatant, add 25 μl 55 mM iodoacetamide to the gel pieces. Allow reaction to proceed in the dark for 45 min at room temperature.
- Remove supernatant. Wash gels with 100 μl NH₄ HCO₃ for 10 min
- Remove supernatant. Wash gels with 100μL of 25 mM NH₄HCO₃ in 50% acetonitrile for 10 min.
- Remove supernatant. Add 100μL of acetonitrile for 10 min to shrink the gels.

- Remove supernatant completely.
- Add trypsin solution to just barely cover the gel pieces. Estimate the gel volume and add about 3x volume of trypsin solution. This volume will vary from sample to sample, but on average 5-25 μ L is sufficient.
- Incubate at 37°C for 4 hours - overnight.
- Supernatant is used for analysis by MALDI TOF MS or ESI QTOF MS.

Contact

Dr. Thomas Ruppert

ZMBH

Im Neuenheimer Feld 282

D-69120 Heidelberg

Phone: +49-6221-54-6895

Fax: +49-6221-54-5891

Email: t.ruppert@zmbh.uni-heidelberg.de