

# GUIDELINES AND INFORMATION FOR CUSTOMERS OF CFMP ZMBH

## Title: Identification of protein interaction partners from pull-down experiments using in-gel approach

Authors: Sabine Merker, Marcin Luzarowski

Reviewed by: Thomas Ruppert

Last updated: 28.07.2023

Contact information: [s.merker@zmbh.uni-heidelberg.de](mailto:s.merker@zmbh.uni-heidelberg.de)

[m.luzarowski@zmbh.uni-heidelberg.de](mailto:m.luzarowski@zmbh.uni-heidelberg.de)

1. PURPOSE
2. DURING INITIAL MEETING INFORM US ABOUT
3. RECOMMENDATIONS
4. PROTEIN ELUTION FROM THE BEADS
5. GENERAL INFORMATION

### 1. PURPOSE

This document describes the recommendations regarding sample preparation and submission to Core Facility for Mass Spectrometry and Proteomics for **identification of interaction partners using pull-down/immunoprecipitation**.

### 2. DURING INITIAL MEETING INFORM US ABOUT

- Do you already have iLab account?  
[https://hmls.corefacilities.org/service\\_center/show\\_external/3564?name=core-facility-for-mass-spectrometry-proteomics](https://hmls.corefacilities.org/service_center/show_external/3564?name=core-facility-for-mass-spectrometry-proteomics)
- Do you work with membrane protein?
- Are you planning to use a detergent during cell lysis, washing or elution of the proteins?
- Did you already perform a test experiment? Did you see proteins on the SDS-PAGE gel when stained with Coomassie?

### 3. RECOMMENDATIONS

- If needed, use only detergents compatible with in-gel digest:  
SDS < 2%  
CHAPS < 4%  
NP-40 < 1%
- You can obtain an aliquot of fixing solution and colloidal Coomassie if you would like to run the gel in your lab but **you can run gel in our facility**. We offer commercial gels (20€ / gel). One can load protein marker and up to **9 samples**. The maximum sample volume is **40 µl**.
- In order to run the gel in our facility, please contact **Sabine Merker by e-mail** and set an appointment.

- Always, bring with you **your elution buffer**. We strongly recommend to load **neighboring lanes** with your **elution buffer**.
- To ensure proper separation on the gel, mix your **marker protein** with appropriate amount of **SDS sample buffer**.
- Use **commercial gels** (higher reproducibility and resolution, reduces contamination)
- For bulk identification of all protein partners, run your gel for **1 cm** (measure the distance from the bottom of the well to the running front). Then simply stop the electrophoresis, wash the gel intensively with dH<sub>2</sub>O and stain with colloidal Coomassie.
- Make a picture of the gel, **mark the lanes on the picture** (preferentially in PowerPoint) and upload it to iLab. File upload is possible after the project request is created.

#### 4. PROTEIN ELUTION FROM THE BEADS

If you work with the antibodies immobilized to the beads (cross-linked), then the most efficient way of releasing the proteins from the beads is to use SDS sample buffer. We recommend using **1x SDS sample buffer**. You can use commercial sample buffer or homemade. Here is the recipe for **4x SDS sample buffer**:

Component	Volume
Tris (1M, pH 6.8)	10 mL
SDS	4 g
B-mercaptoethanol	10 mL
Glycerol	20 mL
Bromophenol blue	0.1 g
dH <sub>2</sub> O	To 50 mL

#### **Elution procedure:**

Mix beads with up to 40 µl of 1x SDS sample buffer. Incubate 10 minutes at 95 °C with shaking. Centrifuge briefly to collect the beads at the bottom of the tube. Transfer the supernatant to the fresh tube.

Note: while working with **membrane proteins**, we recommend to warm sample up to 70 °C or test elution only in presence of SDS (without heating). Membrane proteins tend to aggregate while heated.

#### 5. GENERAL INFORMATION

- We will provide you the results within 4 weeks from the sample submission
- For pull-down experiments, we recommend using at least 3 replicates per condition (e.g. 3 x bait and 3x control).
- Your samples will be analyzed using a 60 min peptide separation method (69€/sample (internal); 86.25€/sample (external)).